

SHAKER FNX 5050



Purpose :

Used to prepare the sample mixture in tubes for Falling Number Analysis.

Features :

- Start and stop buttons
- Mixing time : 3 seconds
- Shakes each tube 20 times
- Stops automatically when test shaking is done.
- Instrument's door has security sensors
- Remove operator influence on results and reduce total analysis time.
- Power: 230V, 50/60 Hz, 6A, 800W.
- Device Dimensions: 220 x 430 x 630 cm, 37 Kg.

Office/ Ofis

Email : export@bastak.com.tr
website: www.bastak.com.tr

Tel : +90 312 395 6 787
Fax: +90 312 395 6 788

İvedik OSB 1354.Cad 1387.Sk No:19
Yenimahalle Ankara TURKEY

Factory / Fabrika

Dağyaka mahallesi 2010 .Cd.No:23
Kazan Ankara TURKEY

Tel : +90 312 811 2 811
Fax: +90 312 811 2 812

Email : info@bastak.com.tr
Website: www.bastak.com.tr



When you work with BASTAK you receive;

- ✓ Free technical consultancy for tailor made solutions,
- ✓ High quality instruments with free webinars and on-line trainings,
- ✓ Fast delivery of instruments, spare parts and consumables,
- ✓ Remote technical support (Audio and Video) six (6) days a week,
- ✓ Economical spare parts and consumables available in 10 years,
- ✓ Free factory visits, live tests and demos,
- ✓ On-Site Warranty and On-Site After Sales Service (for >20.000 USD orders)

“More Than 23 Years of Experience ...”

Bastak Laboratory

57 Types of Device, 35 Types of Additive



Youtube > Bastak Expert

<https://www.youtube.com/watch?v=ktEz2PHdDQg&t=10s>

Bastak Factory - 7000 m2

20.000 instruments in 142 countries



Youtube > Bastak Factory

<https://www.youtube.com/watch?v=UbtWWhEbLy4>



You can find all Catalogues and Competitor Comparison Tables in this link:

https://drive.google.com/drive/folders/1_4aXqljmISG0Gi7c0KeMuPAGWqvi65mT



Sensitive
Nutrition



High
Efficiency



Monoblock
Study Body



Standard
Grinding



Personal
Security

**The information in this document is subject to change without any prior notice or obligation on the part of the manufacturer and should not be construed as a commitment by the manufacturer. The manufacturer assumes no responsibility for any errors that may appear in this document.*

Office/ Ofis

Email : export@bastak.com.tr
website: www.bastak.com.tr

Tel : +90 312 395 6 787
Fax: +90 312 395 6 788

İvedik OSB 1354.Cad 1387.Sk No:19
Yenimahalle Ankara TURKEY

Factory / Fabrika

Dağyaka mahallesi 2010 .Cd.No:23
Kazan Ankara TURKEY

Tel : +90 312 811 2 811
Fax: +90 312 811 2 812

Email : info@bastak.com.tr
Website: www.bastak.com.tr